

Wide-scope analysis of veterinary drug and pesticide residues in animal feed by liquid chromatography coupled to quadrupole-time-of-flight mass spectrometry

María M. Aguilera-Luiz · Roberto Romero-González ·
Patricia Plaza-Bolaños · José Luis Martínez Vidal ·
Antonia Garrido Frenich

Received: 22 February 2013 / Revised: 2 May 2013 / Accepted: 10 May 2013 / Published online: 28 May 2013
© Springer-Verlag Berlin Heidelberg 2013

Abstract A fast and generic method has been developed for the simultaneous monitoring of >250 pesticides and veterinary drugs (VDs) in animal feed. A ‘dilute-and-shoot’ extraction with water and acetonitrile (1 % formic acid) followed by a clean-up step with Florisil cartridges was applied. The extracts were analysed by ultra-high performance liquid chromatography coupled to hybrid analyser quadrupole-time-of-flight mass spectrometry using both positive and negative electrospray ionisation. The detection of the residues was accomplished by retention time and accurate mass using an in-house database. The identification of the detected compounds was carried out by searching of fragment ions for each compound and isotopic pattern. The optimised method was validated and recoveries ranged from 60 % to 120 % at three concentrations (10, 50 and 100 $\mu\text{g kg}^{-1}$) for 30 %, 68 % and 80 % of compounds, respectively, included in the database (364) in chicken feed. Document SANCO 12495/2011 and Directive 2002/657/CE were used as guidelines for method validation. Intra-day and inter-day precisions, expressed as relative standard deviations, were lower than 20 % for more than 90 % of compounds. The limits of quantification ranged from 4 to 200 $\mu\text{g kg}^{-1}$ for most analytes, which are sufficient to verify compliance of products with legal tolerances. The applicability of the

procedure was further tested on different types of feed (chicken, hen, rabbit and horse feed), evaluating recoveries and repeatability. Finally, the method was applied to the analysis of 18 feed samples, detecting some VDs (sulfadiazine, trimethoprim, robenidin and monensin Na) and only one pesticide (chlorpyrifos).

Keywords Animal feed · Pesticides · Veterinary drugs · UHPLC · QqTOF–MS · Dilute-and-shoot

Introduction

The use of pesticides and veterinary drugs (VDs) in food industries is widespread, and they are utilised to fight against plagues [1] and treat diseases that could affect animals [2]. However, their fraudulent use can result in the presence of residues in the final animal products [3, 4], which involves several problems to public health such as the increased risk of allergies, the development of antibiotic-resistant bacteria and harmful effects on the human system (e.g. carcinogenic effects) [5, 6].

The study of pesticide and VD residues in animal feed is an important issue because they can reach the food chain. On the one hand, animal feed is usually employed to administer drugs to animals [7] and, on the other hand, they can be prepared with several ingredients, which are not always free of pesticides [3]. Consequently, in order to ensure animal safety, the European Union (EU) has established regulations for the control of residues in feed, setting maximum residue limits (MRLs) for pesticides [8] and VDs [9, 10]. Moreover, Regulation (EC) 1831/2003 banned the use of all antibiotics other than coccidiostats and histomonostats as feed additives (used as growth

Electronic supplementary material The online version of this article (doi:10.1007/s00216-013-7060-5) contains supplementary material, which is available to authorized users.

M. M. Aguilera-Luiz · R. Romero-González · P. Plaza-Bolaños ·
J. L. Martínez Vidal · A. Garrido Frenich (✉)
Department of Chemistry and Physics (Analytical Chemistry Area),
Research Centre for Agricultural and Food Biotechnology (BITAL),
University of Almería, Agrifood Campus of International
Excellence, ceiA3, 04120 Almería, Spain
e-mail: agarrido@ual.es

promoters) since 1 January 2006. Furthermore, these products have been completely banned on December 2012 [11]. Thus, VD_s should be strictly applied for the prevention or treatment of diseases. In this way, these regulated limits require the determination of an extremely high number of VD and pesticide residues in feed with the aim of ensuring their use in a rational and safe way.

Therefore, sensitive and reliable methodologies have been developed for the analysis of each group of residues separately. The determination of pesticides in feed has been based on Soxhlet [12–14], pressurised liquid extraction (PLE) [15], matrix solid-phase dispersion (MSPD) [16] or QuEChERS [17]. Some of these procedures have been followed by a clean-up step based on solid-phase extraction (SPE) [14, 18] or dispersive SPE (d-SPE) [17], and subsequent chromatographic analysis by gas chromatography (GC) coupled to electron capture detection (ECD) [12, 14], mass spectrometry (MS) [13, 15] or tandem mass spectrometry (MS/MS) [12, 14, 16–18]. For the determination of VD_s, solvent extraction assisted by mechanical shaking [19–22], modified QuEChERS methodology [23] and PLE [24, 25] followed by SPE clean-up [19–21, 24, 25] or liquid–liquid extraction (LLE) [22] have been utilised, using liquid chromatography coupled to MS (LC–MS) [19] or MS/MS (LC–MS/MS) [20–25] methodologies.

Normally, most of these methods are focused on specific groups of residues, not being suitable for wide-scope multi-residue analysis able to cover different classes of residues and contaminants. However, the use of wide-scope methods inherently involves the need for generic extraction procedures, and this strategy shortens the possibilities for clean-up steps [26, 27], in order to minimise analyte losses. Nevertheless, the lack of selectivity in sample preparation can be compensated by selectivity/sensitivity in instrumental analysis. Although in terms of sensitivity, the use of GC or LC coupled to triple quadrupole (QqQ) is usually preferred, and the number of compounds that can be acquired by these techniques in one run is limited [28, 29]. The use of high-resolution mass spectrometry (HRMS), such as time-of-flight (TOF), hybrid quadrupole-time-of-flight (QqTOF) or Orbitrap [29], is considered as a powerful alternative. These full-scan technologies provide the selectivity and the sensitivity required for efficient and wide-range screening, as they combine high, full-spectral sensitivity with high mass resolution, allowing the measurement of accurate mass of any ionisable component in the sample. In this way, the combination of ultra-high performance liquid chromatography (UHPLC) with TOF, QqTOF or Orbitrap analysers has demonstrated to be an excellent tool for the analysis of different classes of residues and contaminants in food and environmental samples [30–32]. However, the number of applications in feed matrices is quite limited [28, 33, 34], where a group of pesticides, VD_s, mycotoxins and plant

toxins are determined in horse feed [28, 33] and coccidiostats in animal feed [34].

Therefore, the aim of this work is the development of a reliable and sensitive methodology based on the use of advanced UHPLC–QqTOF–MS technique for the quantitative determination and reliable identification of more than 250 pesticides and VD_s in different types of feed matrices using generic methods.

Experimental

Reagents and apparatus

VD and pesticide analytical standards were purchased from Riedel-de-Haën (Seelze, Germany), Dr. Ehrenstorfer GmbH (Ausburg, Germany), Sigma-Aldrich (Madrid, Spain), Fluka (Steinheim, Germany), Santa Cruz (Santa Cruz, CA, USA), European Pharmacopoeia (Strasbourg, France), Witega (Berlin, Germany) and LGC Standards (Barcelona, Spain). Stock standard solutions of individual compounds (200–400 mg L^{−1}) were prepared in methanol or acetonitrile (HPLC grade, Fluka), and they were stored at 5 °C or −18 °C (VD_s). A solution for each family of VD_s was prepared from corresponding individual stock standard solutions in methanol or acetonitrile and a multi-pesticide solution was prepared in methanol. Tetracycline and penicillin solutions were renewed monthly. Then, a multi-compound working solution containing all the analytes (0.31 mg L^{−1}) was prepared by combining suitable aliquots of each individual stock standard solution and diluting them with LC–MS-grade methanol, obtained from Fluka. This solution was kept at −18 °C.

Ultrapure water, obtained from a Milli-Q Gradient water system (Millipore, Bedford, MA, USA), was used for the preparation of all aqueous solutions. Formic acid (purity >98 %) and ammonium formate (purity >99 %) were obtained from Panreac (Barcelona, Spain). LC–MS water was purchased from Scharlau (Barcelona, Spain) and LC–MS acetonitrile was obtained from Fluka. ExtraBond Florisil cartridges (500 mg) were provided by Scharlau (Barcelona, Spain). D-SPE sorbents for clean-up experiments included C18-bonded silica, which was purchased from Agilent technologies (Santa Clara, CA, USA), and Florisil from J.T. Backer (Deventer, Holland).

For accurate mass calibration of the QqTOF instrument, a sodium formate solution (prepared from a mixture of 100 µL of 0.1 M sodium hydroxide solution, 200 µL of an aqueous solution of 10 % of formic acid and 20 mL of acetonitrile/water (80:20, v/v)) and a leucine–enkephalin solution at 500 pg µL^{−1} (prepared from 500 µL of 1 mg mL^{−1} leucine–enkephalin solution in water with 1,000 mL of acetonitrile/water 50:50, v/v and 0.1 % formic

acid) were used as calibration and lock-mass solution, respectively. All reagents were purchased from Sigma-Aldrich.

Centrifugations were performed in a high-volume centrifuge Centronic II (JP Selecta, Barcelona, Spain). An analytical AB204-S balance (Mettler Toledo, Greinfensee, Switzerland) was also used. A Reax-2 rotary agitator from Heidolph (Schwabach, Germany) was used for sample extraction.

UHPLC–QqTOF–MS/MS analysis

The analyses were carried out using an Acquity UPLC system using an Acquity UPLC BEH C18 column (2.1 mm×100 mm, 1.7 µm particle size), both from Waters (Milford, MA, USA). The mobile phase consisted of 0.1 % (v/v) formic acid and ammonium formate 4 mM in water (eluent A) and 0.1 % (v/v) formic acid and ammonium formate 4 mM in methanol (eluent B). The analysis started with 95 % of eluent A, which was kept constant during 1 min. Then, it was linearly decreased up to 0 % in 7 min. This composition was held during 4 min before returning to the initial composition in 0.5 min, followed by a re-equilibration time of 1.5 min, giving a total running time of 14 min. The flow rate was 0.3 mL min⁻¹ and the column temperature was set at 30 °C. Aliquots of 5 µL of the sample extract were injected into the chromatographic system.

The UHPLC system was coupled to a QqTOF (XevoTM QtofMS, Waters) equipped with an electrospray interface (ESI) operating in positive (ESI+) and negative ionisation mode (ESI-) using the following parameters: cone voltage, 30 V; capillary voltage, 3.5 kV for ESI (+) and 3.4 kV for ESI (-); desolvation temperature, 350 °C; source temperature, 120 °C; cone gas, 80 L h⁻¹; and desolvation gas, 600 L h⁻¹. Mass spectra were acquired in positive-ion and negative-ion centroid mode. The microchannel plate (MCP) detector was set at 2,000 V. Resolution was at least 9,000 FWHM at *m/z* of the lock mass (see below). Dynamic range enhancement (DRE) was switched on. Mass range in the full scan experiments was set at *m/z* 90–1,000 with a scan time of 0.050 s and interscan time of 0.025 s. In order to perform the fragmentation acquisition for identification purposes, MS/MS (MS^{ETM}) mode was used. This fragmentation mode consists in two simultaneous MS acquisition functions: (1) low energy function (LE) using a collision energy of 4 eV with poor or no fragmentation; and (2) high energy function (HE) using collision energy of 30 eV, promoting collision cell fragmentation. All the analyses were performed using a lock spray with internal lock mass of a solution of leucine-enkephalin (¹²C[M+H]⁺, *m/z* 556.2771 and ¹²C[M-H]⁻, *m/z* 554.2620) delivered to the ESI source at 10 µL min⁻¹. The instrument was verified or calibrated before the analysis by using the calibration solution.

Data acquisition was performed using MassLynx 4.1 software with QuanLynx program (Waters).

Feed samples

Different types of feed samples were purchased from several markets in Almería, Granada, Murcia (Spain) and local farms located in Almería. All samples were grinded and stored at room temperature in the dark. A chicken feed sample was used for optimisation and validation purposes, after checking that it was free of the selected compounds.

Sample extraction

Pesticides and VDs were extracted from feed using an extraction method based on the procedure previously described by Gómez-Pérez et al. [35]. Briefly, 2.5 g of feed was weighed into a 50-mL propylene tube and 2.5 mL of water (LC–MS grade) was added. The mixture was shaken in a vortex. Then 7.5 mL of acetonitrile (1 % formic acid, v/v) was added, and the sample was put into a rotary agitator and extracted by end-over-end shaking for 1 h at 50 rpm. After that, the mixture was centrifuged for 10 min at 4,500 rpm (2,264×*g*), and 1 mL of the supernatant was passed through a Florisil cartridge. The final extract was transferred to a vial and injected into the chromatographic system.

Detection and identification of the analytes

For detection and identification of analytes, an in-house database containing 364 compounds was performed. In order to develop it, a previous characterisation of the analytes was carried out by injection of a multi-compound standard solution containing all the analytes at 250 µg L⁻¹ into the UHPLC–QqTOF–MS system under the conditions previously optimised. Compounds were acquired in both ESI (+) and ESI (-) mode. In this way, after injection in MS mode, the essential information obtained was the retention time (RT), ionisation mode (polarity), characteristic ions and possible adducts. In total, 364 compounds were included in the database to facilitate the analysis of pesticides and VDs in feed matrices. Additional information related to each analyte can be found in the Electronic Supplementary Material (ESM, Table S1).

To identify the compounds, the same multi-compound standard solution was injected into the system employing the MS^{ETM} mode for the detection of the main fragment ions in the HE acquisition function. Thus, one or more fragment ions were selected for each compound and included into the database. In addition to the monitored fragment ions, the identification of the compounds was also carried out by the use of the isotopic pattern.

Validation of the optimised methodology

The proposed method was validated, estimating linearity range, trueness, intra-day and inter-day precision, and limits

of quantification (LOQs). Document SANCO 12495/2011 [36] and Directive 2002/657/CE [37] were used as guidelines for method validation.

Linearity was evaluated performing matrix-matched calibration, using blank feed extracts and adding the corresponding amount of working solution with the target compounds at five concentrations, 1, 2.5, 12.5, 50 and 150 $\mu\text{g kg}^{-1}$. Calibration curves were obtained by least-squares linear regression analysis of the peak area versus concentration. To evaluate trueness, recovery studies were carried out at three concentrations (10, 50 and 100 $\mu\text{g kg}^{-1}$) using spiked blank feed samples ($n=5$), and they were quantified using matrix-matched calibration. Intra-day and inter-day precision (evaluated in five consecutive days) was evaluated at the same concentrations ($n=5$).

The instrument LOQs were estimated by analysing spiked extracted feed samples at low concentrations (from 1 to 50 $\mu\text{g kg}^{-1}$). They were calculated as the minimum concentration for which the accurate mass error was <5 ppm for the characteristic ion of the selected compounds ($[\text{M}+\text{H}]^+$, $[\text{M}+\text{Na}]^+$, $[\text{M}+\text{NH}_4]^+$ or $[\text{M}-\text{H}]^-$) and one fragment ion, as well as the experimental isotopic pattern matched with the theoretical isotopic pattern for these ions.

Results and discussion

In this work, 129 VDs have been studied, including macrolides, quinolones, tetracyclines, sulphonamides, avermectins, nitroimidazoles, coccidiostats, penicillins, amphenicols, tranquilisers, corticoids, ionophores, non-steroidal anti-inflammatory drugs, β -agonists and estrogens, gestagens and androgens. Furthermore, 235 pesticides including organophosphorus, organochlorines, carbamates, neocotinoids, triazines, triazoles and ureas were also selected.

Development of the UHPLC–QqTOF–MS analysis

Firstly, general conditions for the ion source and ion transmission parameters were studied. To this aim, a solvent-based standard at 250 $\mu\text{g L}^{-1}$ was injected into the system, in positive and negative mode, and different parameters of the source and the analyser (scan time, inter-scan time and collision energy of the acquisition functions of LE and HE employed in the MS^E mode) were optimised in order to improve absolute sensitivity of the characteristic ions, observing that the conditions indicated in the previous section provided suitable results.

Secondly, in relation to the chromatographic conditions, generic chromatographic conditions reported in a previous study [35] were used. The overall analysis time was 14 min, and the retention time ranged from 0.9 (cyromazine) to 9.4 min (narsin and salinomycin).

In order to develop a database for the detection of the target compounds, the analytes were characterised and the essential information includes retention time (RT), ionisation mode (polarity), characteristic ions and possible adducts (i.e. with Na^+ or NH_4^+). In this way, a solvent standard solution of the studied compounds at 250 $\mu\text{g L}^{-1}$ was analysed. Because the QqTOF–MS instrument utilised in this study cannot acquire simultaneously in (+) and (–), two separate injections were needed for each polarity mode. The compounds were monitored employing narrow mass extraction windows (10 ppm) within their characteristic RT windows (RTWs). Exact mass values, characteristic ions, RTWs, ionisation mode and mass errors for the studied VDs and pesticides are shown in Table S1, Electronic Supplementary Material.

Regarding the identification of the detected compounds, it was necessary to select fragments for each compound, taking into account the identification requirements provided by the SANCO guidelines [36]. In the case of QqTOF working at high resolution and mass accuracy, the detection of two or more diagnostic ions (quasi-molecular ion and at least one fragment ion) with mass accuracy <5 ppm is required. Besides, the comparison of the obtained and the theoretical isotopic patterns was carried out. For that purpose, a solvent-based standard at 250 $\mu\text{g L}^{-1}$ was injected and MS/MS acquisitions were obtained. In this way, the peak shape of the extracted chromatogram for each fragment (in the HE acquisition function) was compared to that of the main ion (in the LE acquisition function) checking that both traces showed the same RT and elution profile. Moreover, accurate mass errors must be <5 ppm, as well as the isotopic pattern of the fragment ion must be similar to the theoretically profile of the empirical formulae.

Optimisation of the extraction method

The generic extraction method developed by Mol et al. [33] allowed the extraction of a wide set of pesticides, VDs and other contaminants in feed. A slight modification of this method has already been employed for the simultaneous determination of pesticides and VDs in other matrices such as honey in our laboratory [35], and it was selected as the starting point herein. Initially, the experiments were carried out employing blank chicken feed samples spiked at 100 $\mu\text{g kg}^{-1}$. As result of this first assay, 75 % of the compounds (275 compounds) showed acceptable recoveries (acceptable range 60–125 %). Tetracyclines and penicillins were not successfully extracted, and some quinolones and avermectins did not show adequate performance characteristics. This could be due to the high volume of acetonitrile employed, which is not efficient to extract polar analytes like tetracyclines and penicillins because of co-precipitation of some analytes with the proteins. However, the chromatograms obtained after UHPLC–QqTOF–MS analysis showed

numerous background interference peaks for some compounds (Fig. 1a). This can be explained because of the complexity of this matrix, which contains cereals, sugars/fruits, fats, roots/tubers, leguminous/oleaginous plants, amino acids and/or minerals [7]. Therefore, the high complexity of the matrix makes the clean-up of the extracts almost mandatory to obtain reliable results, except in cases where the sensitivity of the method is high enough to offer the option of simply diluting the extracts. Hence, clean-up is an important step when extracts are analysed with TOF, which has showed lower sensitivity than other HRMS analyser, such as Orbitrap, in certain food applications [38]. Bearing in mind this and the results obtained in the first experiment, an additional clean-up stage was checked. In this way, three fast clean-up procedures were evaluated: (1) d-SPE with Florisil, (2) d-SPE with C18 and (3) SPE employing Florisil cartridges. For d-SPE clean-ups, 1 mL of sample extract was taken to an Eppendorf tube containing 100 mg of Florisil or C18 dispersive phase. The mixture was shaken in a vortex for 1 min and then it was centrifuged for 5 min at 4,500 rpm ($2,264\times g$). The final extract was collected and was placed in a vial. To test SPE clean-up, 1 mL of the extract was loaded onto the Florisil cartridge. The final extract was collected and placed in a vial for injection. The obtained results are shown in Fig. 2, and it can be observed that they are similar regarding the number of

extracted compounds, although slightly better recovery results were obtained when SPE-Florisil was used. Furthermore, the clean-up with Florisil cartridges provided 'cleaner' chromatograms (Fig. 1b) with an apparent lower number of interference peaks and slightly better sensibility.

Therefore, the generic approach using a SPE clean-up step was chosen and tested at lower concentrations in the further validation: $10\text{ }\mu\text{g kg}^{-1}$ and $50\text{ }\mu\text{g kg}^{-1}$.

Validation

Ion suppression or matrix enhancement is an important parameter that must be evaluated during method validation. To study matrix effect, two calibration sets (from 1 to $150\text{ }\mu\text{g kg}^{-1}$) were compared: solvent and matrix-matched standard calibration, using chicken feed. Matrix effect was calculated as described elsewhere [35]. It can be noticed that matrix effect was observed for 86 % of the compounds: 49 % of the analytes showed ion suppression and 37 % of the compounds matrix enhancement. On the other hand, only 14 % of the compounds showed a slight ion suppression (7 %) or ion enhancement (7 %). Consequently, matrix-matched standard calibration was used for quantification purposes.

Linearity was then evaluated in the range from 1 to $150\text{ }\mu\text{g kg}^{-1}$. Good linearity was obtained when peak area




Fig. 1 UHPLC–QqTOF–MS chromatograms of both quantification and detection transition of albendazole and meclofenamic acid obtained: **a** without clean-up; **b** clean-up by SPE using Florisil cartridges




Fig. 2 Recovery results obtained for the extraction of spiked chicken feed samples ($100 \mu\text{g kg}^{-1}$) without clean-up step and using different procedures of clean-up, d-SPE procedure employing Florilil, C18 and SPE procedure employing Florilil cartridges

was used as the analytical response, and determination coefficients (R^2) were higher than 0.9800.

Trueness was estimated through recovery studies at three concentrations (10, 50 and $100 \mu\text{g kg}^{-1}$). Figure S1 in the Electronic Supplementary Material shows the total ion chromatograms obtained for the four different feed samples (chicken, hen, rabbit and horse) spiked at $50 \mu\text{g kg}^{-1}$. Table 1 shows a summary of the obtained results, whereas full results can be found in Table S2, Electronic Supplementary Material. It can be observed that 30 % of compounds (111 compounds) showed adequate recovery (recovery ranging from 60 % to 125 %) at the lowest concentration checked ($10 \mu\text{g kg}^{-1}$). As expected, the number of compounds showing acceptable recovery at $50 \mu\text{g kg}^{-1}$ and $100 \mu\text{g kg}^{-1}$ was considerably higher: 68 % and 80 % (247 and 291 compounds), respectively. Intra-day and inter-day precision values (expressed as relative standard derivation, RSD, %) obtained for each concentration are shown

in Table S2, Electronic Supplementary Material. It can be observed that intra- and inter-day precision values were always ≤ 20 % for the majority of compounds (Table S2, Electronic Supplementary Material).

LOQs obtained for the studied compounds (Table S2, Electronic Supplementary Material) were lower than or equal to $10 \mu\text{g kg}^{-1}$ for 43 % of the compounds, equal to $50 \mu\text{g kg}^{-1}$ for 37 % and equal to or higher than $200 \mu\text{g kg}^{-1}$ in sample for 6 % of the target compounds (Table 1). Some compounds showed LOQs higher than the lowest level of concentration at which they can be recovered, as heptenophos or isoxicam (Table S2, Electronic Supplementary Material). It can be explained due to the quasi-molecular ion, which was used for quantification, could be perfectly detected at concentrations lower than the LOQ. However, the fragments were not detected at such low concentration because of the low sensitivity of the fragments. In these cases, the isotopic pattern was used for the reliable identification of the compounds at concentrations lower than LOQ.

Additionally, the feasibility of the method for the analysis of target compounds in different types of feed was checked. Thus, three additional types of feed (hen, rabbit and horse feed) were submitted to validation. Firstly, we tested if chicken feed could be selected as representative matrix for these three types of feed. Matrix effect was evaluated and the same concentrations (from 1 to $150 \mu\text{g kg}^{-1}$) were analysed in each type of feed. Then, the slopes of the calibration curves were compared with the slope of the calibration curve obtained in chicken feed samples. A significant difference was observed among the tested feed matrices. In the case of hen feed, 21 % of the compounds showed significant differences with chicken feed, whereas in rabbit and horse feed, the number of compounds with significantly different slopes was higher, 69 % and 67 %, respectively. Therefore, it was not possible to select one type of feed as a representative matrix, which is probably due to the different composition of each type of feed (type of nutritional properties), requiring specific matrix-matched calibration standards. Moreover, trueness and repeatability were also evaluated following the same procedure employed

Table 1 Summary of the validation results showing percentage of analytes with acceptable recoveries (60–125 %) and LOQs

Feed	% Compounds with acceptable recovery			% Compounds ^a	
	$10 \mu\text{g kg}^{-1}$	$50 \mu\text{g kg}^{-1}$	$100 \mu\text{g kg}^{-1}$	$\text{LOQ} \leq 10 \mu\text{g kg}^{-1}$	$\text{LOQ} > 10 \mu\text{g kg}^{-1}$
Chicken	30 ^a	68	80	43	43 ^b
Hen	25	53	61	28	59
Rabbit	22	51	59	16	71
Horse	24	57	64	17	70

^a Percentage results based on all target compounds initially included in the study (364 compounds)

^b The remaining percentage corresponds with the number of compounds not extracted

in chicken feed, showing the obtained results in Table S3, Electronic Supplementary Material, and a summary in Table 1. It can be observed that for hen feed, the percentage of compounds with adequate recovery and repeatability values was 25 % (90 compounds) at $10 \mu\text{g kg}^{-1}$, 53 % (193 compounds) at $50 \mu\text{g kg}^{-1}$ and 61 % (221 compounds) at $100 \mu\text{g kg}^{-1}$. For rabbit feed, these percentages are 22 % (81 compounds) at $10 \mu\text{g kg}^{-1}$, 51 % (187 compounds) at $50 \mu\text{g kg}^{-1}$ and 59 % (216 compounds) at $100 \mu\text{g kg}^{-1}$. Finally for horse feed, the number of compounds with adequate figures of merit was 24 % (88 compounds) at $10 \mu\text{g kg}^{-1}$, 57 % (209 compounds) at $50 \mu\text{g kg}^{-1}$ and 64 % (234 compounds) at $100 \mu\text{g kg}^{-1}$. As expected, the number of compounds extracted in hen, rabbit and horse feed with acceptable recovery was considerably lower than in chicken feed due to the fact that these matrices have different composition (especially different fatty content) and they can be considered more complex.

Finally, LOQs obtained for each compound in the matrices evaluated are shown in Table S3, Electronic Supplementary Material. It can be observed that 69 % of compounds showed LOQs $\leq 50 \mu\text{g kg}^{-1}$ in hen feed, whereas in rabbit and horse feed, the number of compounds with LOQs $\leq 50 \mu\text{g kg}^{-1}$ was 55 % and 61 %, respectively. Table 1 shows a summary of the results obtained for all the matrices. As aforementioned, some compounds showed LOQs higher than the lowest concentration at which they can be recovered. Again, these compounds showed problems to detect the fragment at LOQ concentration, although their quasi-molecular ion was perfectly detected.

Despite the higher levels of LOQ obtained for some compounds, many of them were lower than MRLs established by the legislation for some specific compounds, as it can be seen in Table 2.

In general, it can be observed that the obtained LOQs were relatively high. This can be explained bearing in mind the acquisition mode of the fragments based on the MS^E procedure. Despite its aforementioned advantages, the success of the MS^E approach can be limited by the quality of the spectrum [31]. The low sensitivity or strongly interfered spectra make difficult the identification of compounds using fragments at low concentrations.

As aforementioned, many analytical methods have been developed for the determination of compounds belonging to one family of pesticides or VDs in feed [7, 14, 23, 39, 40]. Nevertheless, the number of multi-class or wide-scope multi-class methods is scarce. In consequence, it is important to highlight the ability of the proposed method to analyse a high variety of pesticides and VDs in different feed samples.

If the proposed method is compared with the published multi-class methods developed for the determination of VDs [22, 24, 41] and pesticides [15, 17] in feed, it can be indicated that the combination of a generic extraction method with UHPLC–QqTOF–MS analysis is simpler and faster than

Table 2 Comparison of MRL established for some VDs and pesticides in animal feed with LOQs obtained with the proposed method

Compound	Animal feed	MRL ($\mu\text{g kg}^{-1}$)	LOQ ($\mu\text{g kg}^{-1}$)
Diclazuril	Chicken	10	12.5
	Horse	30	50
Endosulfan sulphate	Chicken	100	12.5
Monesin Na	Chicken	1,250	12.5
	Rabbit	1,250	2.5
	Horse	1,250	2.5
Narasin	Hen	700	12.5
	Rabbit	700	12.5
	Horse	2,100	12.5
Robenidine	Hen	700	12.5
	Rabbit	6,600	12.5
	Horse	2,100	12.5
Salinomycin	Chicken	700	12.5
	Hen	700	150
	Rabbit	700	150
	Horse	2,100	50

previous methods. Moreover, the presented methodology increased sample throughput because it only requires a rapid clean-up step through Florisil cartridges instead of d-SPE clean-up [17, 41], LLE followed by concentration step [22] or an extraction by SPE [24]. In relation to the number of validated compounds, the proposed method allowed the determination and quantification of an average of 200 compounds in different feed improving the number of compounds of other works: 33 VDs [41], 10 prohibited medical additives (three analytes were analysed in qualitative mode) [22], 18 antibacterials [24] and 33 [15] and 122 [17] pesticides.

Concerning the performance characteristics, the developed method allowed the extraction and quantification of a high number of VDs and pesticides at low levels ($10 \mu\text{g kg}^{-1}$). These results are better than those obtained in a previous multi-class method developed for the determination of 122 pesticides [17], where some pesticides were validated at a higher concentration, $50 \mu\text{g kg}^{-1}$ (for example atrazine, carbaryl or metalaxyl). In another study, $100 \mu\text{g kg}^{-1}$ was the lowest level at which antibacterial compounds were determined [22]. In contrast to this, a multi-class method developed for the analysis of 33 VDs [41] showed good recoveries at low levels ($<10 \mu\text{g kg}^{-1}$).

Regarding LOQs, the obtained values in the proposed method (ranged from 4 to $200 \mu\text{g kg}^{-1}$ in sample) were lower than those obtained in two previous works [17, 22]. The first one determined pesticides in wheat grain ($7\text{--}400 \mu\text{g kg}^{-1}$) [17] and the second one measured VDs in pig and poultry compound feed ($100\text{--}169 \mu\text{g kg}^{-1}$) [22]. However, other works showed better LOQs, such as those

obtained by Boscher et al. [41] (LOQs ranging from 3.8 to 65 $\mu\text{g kg}^{-1}$ in sample, piglet feed) and by Kantiani et al. [24] (instrumental LOQs from 0.57 to 2.92 $\mu\text{g kg}^{-1}$). These low LOQ values can be explained due to the fact that both works employed QqQ as analyser, which showed better sensitivity than QqTOF using MS^E. Moreover, these methods were developed for a limited number of compounds.

One of the main disadvantages of the proposed method is that it was not able to extract and quantify some families of compounds, such as tetracyclines, quinolones or penicillins, which were analysed in other works [24, 41] with good recoveries, and therefore, current work is focused on this.

If it is compared with the two wide-scope multi-class methods previously developed for analysis of VDs and pesticides in feed samples (horse feed and cereal-based product, respectively) [33, 42], some advantages can be mentioned. Firstly, the generic extraction method was developed from Mol et al. [33], which is similar to the method proposed in this work, with a slight modification by the addition of a simple clean-up step. Secondly, the generic extraction method is faster than the method developed by van der Lee et al. [42], where the extraction time was 18 h and it required a clean-up procedure by gel permeation chromatography and SPE. Furthermore, the number of pesticides and VDs validated with our method was higher or similar than the number of compounds of both works, 106 compounds [42] and 222 [33], respectively. Finally, regarding LOQs, both methods showed similar values compared with our method. In one of these works, the instrumental LOQs ranged from 1 to 46 $\mu\text{g kg}^{-1}$ [42], and in the other one, LODs ranged from 10 to 50 $\mu\text{g kg}^{-1}$ [33].

Analysis of real samples

The optimised method was applied to 18 samples of feed, including rabbit (four), chicken (seven), hen (two), pig (one), lamb (one) and horse (one) feed. Moreover, one sample of

barley and one sample of maize (two components of animal feed) were also selected.

An internal quality control was carried out for every batch of samples to check if the system is under control, and it implies four matrix-matched calibrations (in chicken, hen, rabbit and horse), reagent blanks, full procedural blanks and fortified extracts at 150 $\mu\text{g kg}^{-1}$. Instrument blanks were composed of methanol and were analysed at the beginning of the batch. Since analyte signals were not found in the instrumental blanks, no further actions were taken. Full procedural blanks were prepared by four different types of feed sample extracted and analysed after instrumental blanks. This experiment showed zero carry-over for all the analytes, ensuring that no cross-contamination was taking place.

As aforementioned, analyte detection and identification was based on the criteria showed in section ‘[Detection and identification of the analytes](#)’.

When the samples were analysed by the proposed method, seven samples contained traces of VDs and pesticides (Table 3). Some VDs (sulfadiazine, trimethoprim, robenidine and monensin sodium) and one pesticide (chlorpyrifos) were detected. Sulfadiazine and trimethoprim were detected in two chicken samples. One of them contained 1,114 $\mu\text{g kg}^{-1}$ of sulfadiazine and 225 $\mu\text{g kg}^{-1}$ of trimethoprim, and the other one contained 217 $\mu\text{g kg}^{-1}$ of sulfadiazine and 72 $\mu\text{g kg}^{-1}$ of trimethoprim. Sulfadiazine, which is an antibiotic that belongs to the group of sulphonamides, is a synthetic bacteriostatic compound that is commonly used to treat urinary tract infections. On the other hand, trimethoprim is a dihydropteroate synthetase inhibitor that is commonly used in combination with sulphonamides for broad-spectrum antimicrobial therapy. In veterinary medicines, the combination of sulphonamides with trimethoprim is widely used for prevention and treatment of respiratory or gastro-intestinal tract infections in cattle, swine and poultry [43]. Robenidine was detected in a rabbit sample at 4,186 $\mu\text{g kg}^{-1}$, being this level lower than the MRL

Table 3 Concentration of VDs and pesticides in the analysed samples

Compound	Concentration ($\mu\text{g kg}^{-1}$)							MRL ^c
	M1 ^a	M2	M5	M11	M15	M16	M17	
Chlorpyrifos	65^b	18				92	193	10
Sulfadiazine			1,114	217				N.E. ^d
Trimethoprim			225	72				N.E.
Robenidin		4,186						6,600
Monensin Na	124	315			239	189	84	1,250


^a Type of feed is indicated herein: rabbit, M1 and M2; chicken, M5 and M11; pig, M15; lamb, M16; horse, M17

^b Figures in bold indicate that the concentration in sample exceeded MRL

^c MRL maximum residue level

^d N.E. not established

Fig. 3 UHPLC–QqTOF–MS/MS extracted chromatograms of a positive sample of rabbit feed showing quantification and detection transition (*upper*) and confirmation (*lower*) transitions containing **a** monensin sodium at $315 \mu\text{g kg}^{-1}$ and **b** robenidine at $4,186 \mu\text{g kg}^{-1}$



established in legislation ($6,600 \mu\text{g kg}^{-1}$) [44]. Its presence in this sample can be explained due to this coccidiostat can be used in fattening and breeding rabbits in feedstuffs [45]. The coccidiostat monensin sodium was detected in five feed samples: rabbit (two), pig (one), lamb (one) and horse (one) feed (Table 3). This coccidiostat is used as additive in feed to protect animals against coccidiosis in animal husbandry, which can explain its occurrence. However, all concentrations (up to $315 \mu\text{g kg}^{-1}$) were below the established MRLs for each matrix. Figure 3 shows the extracted ion chromatograms for sample 2 (rabbit feed sample) containing $4,186 \mu\text{g kg}^{-1}$ of robenidine and $315 \mu\text{g kg}^{-1}$ of monensin sodium.

An important result concerns that only one pesticide was found, chlorpyrifos, which was detected in four samples at concentrations higher than the MRL (Table 3) in all cases. Its occurrence can be explained as a consequence of its wide application in agriculture and, thus, in the cereals used to produce feed. Besides, this pesticide is defined by the United States Environmental Protection Agency (EPA) as ‘one of the most widely used organophosphate insecticides’ [46].

Conclusions

As aforementioned, the development of wide-scope multi-residue methods able to cover the analysis of pesticide and VD residues in feed is a current need due to the gap that exists in this area of feed safety. Consequently, in this work, an UHPLC–QqTOF–MS method has been developed for

simultaneous determination of a huge number of target pesticides and VD residues in feed. A generic sample procedure followed by a simple clean-up step by Florisil cartridges was selected as the most appropriate for the analysis of the wide range of compounds in feed in terms of recovery yields and matrix effect. The developed method allows the analysis of an average of 200 residues of different families of pesticides and VDs in various feed. The developed method was validated according to Directive 2002/657/CE and Document SANCO 12495/2011; it showed good quantitative results for most of the studied compounds and the LOQs range from 4 to $200 \mu\text{g kg}^{-1}$. Despite the variety of compounds that could be determined by the current method, further experiments are needed to include the analysis of tetracyclines, quinolones, avermectins and penicillins, which could not be determined with this method. Although the sensitivity attained by the proposed method was lower than some of published methods, it was enough to comply with MRLs established by EU. Finally, the method was applied to different kinds of feed demonstrating its applicability as quantitative method in any type of feed. This method illustrated the feasibility of using QqTOF for detecting, identifying and quantifying pesticides and VDs residues in a multi-class method. In spite of the advantages showed by the proposed method, it shows some limitations due to some families of VDs, such as tetracyclines, quinolones and penicillins, which are not properly extracted. Furthermore, the proposed method showed low sensitivity making the detection of some compounds and fragments at low concentrations unfeasible.

Acknowledgements The authors gratefully acknowledge the Spanish Ministry of Economy and Competitiveness (MINECO) and FEDER for financial support (Project Ref. AGL2010-21370). M.M.A.L. acknowledges her grant (F.P.U.) (MINECO, Ref. AP2008-02811). R.R.G. is also grateful for personal funding through the Ramón y Cajal Program (MINECO-ESF). P.P.B. gratefully acknowledges personal funding through the Agrifood Campus of International Excellence, ceiA3 (Spanish Ministry of Education, Culture and Sport).

References

- Wang J, Chow W, Cheung W (2011) Application of a tandem mass spectrometer and core-shell particle column for the determination of 151 pesticides in grains. *J Agric Food Chem* 59:8589–8608
- Stolker AAM, Brinkman UAT (2005) Analytical strategies for residue analysis of veterinary drugs and growth-promoting agents in food-producing animals—a review. *J Chromatogr A* 1067:15–53
- MacLachlan DJ, Bhula R (2008) Estimating the residue transfer of pesticides in animal feedstuffs to livestock tissues, milk and eggs: a review. *Aust J Exp Agric* 48:589–598
- Garrido Frenich A, Plaza-Bolaños P, Aguilera-Luiz MM, Martínez Vidal JL (2010) In: Quintin TJ (ed) *Chromatography types, techniques and methods*. Nova, New York
- Moreno-Bondi MC, Marazuela MD, Herranz S, Rodriguez E (2009) An overview of sample preparation procedures for LC–MS multiclass antibiotic determination in environmental and food samples. *Anal Bioanal Chem* 395:921–946
- Ritter L, Gousheff NCI, Arbuckle T, Cole D, Raizenne M (2006) Addressing the linkage between exposure to pesticides and human health effects—research trends and priorities for research. *J Toxicol Environ Health B Crit Rev* 9:441–456
- Borrás S, Companyó R, Granados M, Guiteras J, Pérez-Vendrell AM, Brufau J, Medina M, Bosch J (2011) Analysis of antimicrobial agents in animal feed. *Trends Anal Chem* 30:1250–1268
- Regulation (EC) No 299/2008 of the European Parliament and of the Council of 11 March 2008 amending Regulation (EC) No 396/2005 on maximum residue levels of pesticides in or on food and feed of plant and animal origin, as regards the implementing powers conferred on the Commission, Official Journal of the European Union L 97/67; available at <http://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2008:097:0067:0071:en:PDF>. Accessed 25 Apr 2013
- Commission Directive 2009/8/EC of 10 February 2009 amending Annex I to Directive 2002/32/EC of the European Parliament and of the Council as regards maximum levels of unavoidable carry-over of coccidiostats or histomonostats in nontarget feed, Official Journal of the European Union L 40/19; available at <http://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2009:040:0019:0025:EN:PDF>. Accessed 25 Apr 2013
- Directive 2002/32/EC of the European Parliament and of the Council of 7 May 2002 on undesirable substances in animal feed, L 140/10. (consolidated version: 2002 L0032—ES—01.01.2012—014.002—1); available at <http://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2002:140:0010:0021:ES:PDF>. Accessed 25 Apr 2013
- Regulation (EC) No 1831/2003 of the European Parliament and of the Council of 22 September 2003 on additives for use in animal nutrition, Official Journal of the European Union L 268/29; available at <http://imm.jrc.ec.europa.eu/CollectionDocuments/EC-1831-2003.pdf>. Accessed 25 Apr 2013
- Aulakh RS, Gill JPS, Bedi JS, Sharma JK, Joia BS, Ockerman HW (2006) Organochlorine pesticide residues in poultry feed, chicken muscle and eggs at a poultry farm in Punjab, India. *J Sci Food Agric* 86:741–744
- Guo Y, Yu HY, Zhang BZ, Zeng EY (2009) Persistent halogenated hydrocarbons in fish feeds manufactured in South China. *J Agric Food Chem* 57:3674–3680
- Iglesias-García I, Barriada-Pereira M, González-Castro MJ, Muniategui-Lorenzo S, López-Mahía P, Prada-Rodríguez D (2008) Development of an analytical method based on microwave-assisted extraction and solid phase extraction cleanup for the determination of organochlorine pesticides in animal feed. *Anal Bioanal Chem* 391:745–752
- Fernandez-Alvarez M, Llompart M, Lamas JP, Lores M, García-Jares C, García-Chao M, Dagnac T (2009) Simultaneous extraction and cleanup method based on pressurized solvent extraction for multiresidue analysis of pesticides in complex feed samples. *J Agric Food Chem* 57:3963–3973
- Carro AM, Lorenzo RA, Fernández F, Rodil R, Cela R (2005) Multi-residue screening of chlorinated and brominated compounds from aquaculture samples using matrix solid-phase dispersion–gas chromatography–mass spectrometry. *J Chromatogr A* 1071:93–98
- Walorczyk S (2007) Development of a multi-residue screening method for the determination of pesticides in cereals and dry animal feed using gas chromatography–triple quadrupole tandem mass spectrometry. *J Chromatogr A* 1165:200–212
- Nácher-Mestre J, Serrano R, Portolés-Nicolau T, Hernández F, Benedito-Palos L, Pérez-Sánchez J (2009) A reliable analytical approach based on gas chromatography coupled to triple quadrupole and time-of-flight mass analyzers for the determination and confirmation of polycyclic aromatic hydrocarbons in complex matrices from aquaculture activities. *Rapid Commun Mass Spectrom* 23:2075–2086
- Van Poucke C, De Keyser K, Baltusnikienė A, McEvoy JDG, Van Peteghem C (2003) Liquid chromatographic–tandem mass spectrometric detection of banned antibacterial growth promoters in animal feed. *Anal Chim Acta* 483:99–109
- Vincent U, Chedin M, Yasar S, von Holst C (2008) Determination of ionophore coccidiostats in feeding stuffs by liquid chromatography–tandem mass spectrometry. Part I. Application to targeted feed. *J Pharm Biomed Anal* 47:750–757
- van Holthoon F, Mulder PJP, van Bennekom EO, Heskamp H, Zuidema T, van Rhijn HJA (2010) Quantitative analysis of penicillins in porcine tissues, milk and animal feed using derivatisation with piperidine and stable isotope dilution liquid chromatography tandem mass spectrometry. *Anal Bioanal Chem* 396:3027–3040
- Cronly M, Behan P, Foley B, Malone E, Earley S, Gallagher M, Shearan P, Regan L (2010) Development and validation of a rapid multi-class method for the confirmation of fourteen prohibited medicinal additives in pig and poultry compound feed by liquid chromatography–tandem mass spectrometry. *J Pharm Biomed Anal* 53:929–938
- Cronly M, Behan P, Foley B, Malone E, Shearan P, Regan L (2011) Determination of eleven coccidiostats in animal feed by liquid chromatography–tandem mass spectrometry at cross contamination levels. *Anal Chim Acta* 700:26–33
- Kantiani L, Farré M, Grases IFJM, Barceló D (2010) Development and validation of a pressurised liquid extraction liquid chromatography–electrospray–tandem mass spectrometry method for β -lactams and sulfonamides in animal feed. *J Chromatogr A* 1217:4247–4254
- Kantiani L, Farré M, Grases IFJM, Barceló D (2010) Determination of antibacterials in animal feed by pressurized liquid extraction followed by online purification and liquid chromatography–electrospray tandem mass spectrometry. *Anal Bioanal Chem* 398:1195–1205
- Stolker AAM, Rutgers P, Oosterink E, Lasaroms JJP, Peters RJB, van Rhijn JA, Nielen MWF (2008) Comprehensive screening and quantification of veterinary drugs in milk using UPLC–ToF–MS. *Anal Bioanal Chem* 391:2309–2322

27. Kinsella B, O'Mahony J, Malone E, Moloney M, Cantwell H, Furey A, Danahera M (2009) Current trends in sample preparation for growth promoter and veterinary drug residue analysis. *J Chromatogr A* 1216:7977–8015
28. Kellmann M, Muenster H, Zomer P, Mol H (2009) Full scan MS in comprehensive qualitative and quantitative residue analysis in food and feed matrices: how much resolving power is required? *J Am Soc Mass Spectrom* 20:1464–1476
29. Kaufmann A (2012) The current role of high-resolution mass spectrometry in food analysis. *Anal Bioanal Chem* 403:1233–1249
30. Ortelli D, Cognard E, Jan P, Edler P (2009) Comprehensive fast multiresidue screening of 150 veterinary drugs in milk by ultra-performance liquid chromatography coupled to time of flight mass spectrometry. *J Chromatogr B* 877:2363–2374
31. Diaz R, Ibáñez M, Sancho JV, Hernández F (2012) Target and non-target screening strategies for organic contaminants, residues and illicit substances in food, environmental and human biological samples by UHPLC–QTOF–MS. *Anal Methods* 4:196–209
32. Nurmi J, Pellinen J (2011) Multiresidue method for the analysis of emerging contaminants in wastewater by ultra performance liquid chromatography–time-of-flight mass spectrometry. *J Chromatogr A* 1218:6712–6719
33. Mol HGJ, Plaza-Bolaños P, Zomer P, de Rijk TC, Stolker AAM, Mulder PPJ (2008) Toward a generic extraction method for simultaneous determination of pesticides, mycotoxins, plant toxins, and veterinary drugs in feed and food matrices. *Anal Chem* 80:9450–9459
34. van der Heeft E, Bolck YJC, Beumer B, Nijrolder AWJM, Stolker AAM, Nielen MWF (2009) Full-scan accurate mass selectivity of ultra-performance liquid chromatography combined with time-of-flight and orbitrap mass spectrometry in hormone and veterinary drug residue analysis. *J Am Soc Mass Spectrom* 20:451–463
35. Gómez-Pérez ML, Plaza-Bolaños P, Romero-González R, Martínez-Vidal JL, Garrido-Frenich A (2012) Comprehensive qualitative and quantitative determination of pesticides and veterinary drugs in honey using liquid chromatography–orbitrap high resolution mass spectrometry. *J Chromatogr A* 1248:130–138
36. Document No SANCO/12495/2011. Method validation and quality control procedures for pesticide residues analysis in food and feed
37. Commission Decision of 12 August 2002 implementing Council Directive 96/23/EC concerning the performance of analytical methods and the interpretation of results, Official Journal of the European Union L 221/8; available at http://eurl.craw.eu/img/page/Legislation/657_2002_EN.pdf Accessed 25 Apr 2013
38. Peters RJB, Stolker AAM, Mol JGJ, Lommen A, Lyris E, Angelis Y, Vonaparti A, Stamou M, Georgakopoulos C, Nielen MWF (2010) Screening in veterinary drug analysis and sports doping control based on full-scan, accurate-mass spectrometry. *Trends Anal Chem* 29:1250–1268
39. van Holthoorn F, Mulder PPJ, van Bennekom EO, Heskamp H, Zuidema T, van Rhijn HJA (2010) Quantitative analysis of penicillins in porcine tissues, milk and animal feed using derivatisation with piperidine and stable isotope dilution liquid chromatography tandem mass spectrometry. *Anal Bioanal Chem* 396:3027–3040
40. Juan C, Igualada C, Moragues F, León N, Mañes J (2010) Development and validation of a liquid chromatography tandem mass spectrometry method for the analysis of β -agonists in animal feed and drinking water. *J Chromatogr A* 1217:6061–6068
41. Boscher A, Guignard C, Pellet T, Hoffmann L, Bohn T (2010) Development of a multi-class method for the quantification of veterinary drug residues in feedingstuffs by liquid chromatography–tandem mass spectrometry. *J Chromatogr A* 1217:6394–6404
42. van der Lee MK, van der Weg G, Traag WA, Mol HGJ (2008) Qualitative screening and quantitative determination of pesticides and contaminants in animal feed using comprehensive two-dimensional gas chromatography with time-of-flight mass spectrometry. *J Chromatogr A* 1186:325–339
43. de Paula FCCR, de Pietro AC, Cass QB (2008) Simultaneous quantification of sulfamethoxazole and trimethoprim in whole egg samples by column-switching high-performance liquid chromatography using restricted access media column for on-line sample clean-up. *J Chromatogr A* 1189:221–226
44. Commission implementing regulation (EU) 532/2011 of 31 May 2011 concerning the authorisation of robenidine hydrochloride as a feed additive for rabbits for breeding and rabbits for fattening (holder of authorisation Alpharma Belgium BVBA) and amending Regulations (EC) No 2430/1999 and (EC) No 1800/2004, Official Journal of the European Union L 146/7; available at <http://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2011:146:0007:0010:EN:PDF>. Accessed 25 Apr 2013
45. Vancraeynest D, De Gussem M, Marien M, Maertens L. The anticoccidial efficacy of robenidine hydrochloride in *Eimeria* challenged rabbits. Pathology and Hygiene, 9th World Rabbit Congress—June 10–13, 2008—Verona—Italy, 1103–1106
46. Interim Reregistration Eligibility Decision for Chlorpyrifos, United States Environmental Protection Agency, Prevention, Pesticides and Toxic substances, EPA 738-R-01-007; available at http://www.epa.gov/oppsrrd1/REDs/chlorpyrifos_ired.pdf. Accessed 25 Apr 2013